

The BlueLine Electrode Product Line

BlueLine 11pH	Combination pH	Silamid	Platinum	Cone	3M KCl	0-14	0-100	No	Standard	Plug Head
BlueLine 12pH	Combination pH	Silamid	Platinum	Cone	3M KCl	0-14	0-100	No	Standard	DIN
BlueLine 13pH	Combination pH	Silamid	Ground Joint	Sphere	3M KCl	0-14	0-100	No	Slurries	Plug Head
BlueLine 14pH	Combination pH	Silamid	Platinum	Cone	3M KCl	0-14	0-100	NTC	3-in-2	DIN + Banana
BlueLine 15pH	Combination pH	Silamid	Platinum	Cone	3M KCl	0-14	0-100	NTC	3-in-2	BNC + Banana
BlueLine 16pH	Combination pH	Silamid	Platinum	Cylinder	3M KCl	0-14	0-100	No	Micro	Plug Head
BlueLine 17pH	Combination pH	Silamid	Platinum	Cone	3M KCl	0-14	0-100	No	Standard	BNC
BlueLine 18pH	Combination pH	Silamid	Platinum	Cone	3M KCL	0-14	(-)5-100	PT1000	3-in-2	DIN + Banana
BlueLine 21pH	Combination pH	Silamid	Hole	Spear	Referid	2-13.	(-)5-80	No	Stabbing	Plug Head
BlueLine 22pH	Combination pH	Silamid	Fibre	Cylinder	Gel	0-14	(-)5-80	No	Noryl	Plug Head
BlueLine 23pH	Combination pH	Silamid	Fibre	Cylinder	Gel	0-14	(-)5-80	No	Noryl	DIN
BlueLine 24pH	Combination pH	Silamid	Fibre	Cylinder	Gel	0-14	(-)5-80	NTC	3-in-2/Noryl	DIN + Banana
BlueLine 25pH	Combination pH	Silamid	Fibre	Cylinder	Gel	0-14	(-)5-80	No	Noryl	BNC
BlueLine 26pH	Combination pH	Silamid	Fibre	Cylinder	Gel	0-14	(-)5-80	NTC	3-in-2/Noryl	BNC + Banana
BlueLine 27pH	Combination pH	Silamid	Annular Gap	Flat	Referid	2-13.	(-)5-50	No	Flat Surfaces	Plug Head
BlueLine 31Rx	Combination ORP	Silamid	Ceramic	Disk	3M KCl	N/A	0-100	No	Redox	Plug Head